

SECTION 3

EDUCATION

Learning is a lifelong endeavour. From simulation products upon which to practice and perfect new techniques to continuing education workshop and lectures that challenge and inform, Gunz Dental can assist.

For all your Education needs Talk to Us.



SECTION 3

EDUCATION

CONTENTS

FRASACO PROFESSIONAL EDUCATION SIMULATION PRODUCTS	39
---	----

FRASACO PROFESSIONAL EDUCATION



Simulation Products

FRASACO

frasaco offers the most comprehensive range of dental simulators, study models and training models for hands-on practice in pre-clinical and post-graduate students' laboratories in preventive and operative dentistry as well as fixed, removable and full prosthodontics, paediatric dentistry and orthodontics, radiology and anesthesia, endodontics, exodontia and oral surgery, oral hygiene, periodontology, gerodontology and implantology, pathology, anatomy and dental technique.

Here is an overview of the most popular frasaco simulation products. For a complete view of the frasaco range, technical specifications or to discuss custom projects such as special teeth or models, please contact your Gunz representative or request a frasaco catalogue.



Simulation Products

Introducing the new ANA-4 and ANA-4V models:

For decades, frasaco dental models have been the first choice in leading dental schools worldwide. Now, frasaco presents the all new ANA-4 Standard model, sure to be the perfect tool for tomorrow's dental education. Designed to represent the typical natural adult dentition, this advanced model was developed from an actual patient impression and was further refined with input and feedback from many instructors. Utilising the latest CAD/CAM technology from design to production, the new-generation ANA-4 (Fig 1.) ensures quality, precision, and repeatability. Please contact us for more information on this 21st century model.

Known throughout the world for its detail, versatility, and realism, the Frasaco AG-3 (Fig 2.) provides an excellent platform for teaching and learning fundamentals of dentistry.

Widely used in the study of Periodontics, as well as Dental Hygiene, the A-PZ (Fig 3.) features realistic examples of periodontitis in the early to advanced stages.

The first and foremost modular-style phantom head on the market, the P-6/3 (Fig 4.) has been a primary tool for teaching dental students procedural skills and techniques for over 20 years. Simulated jaw fully articulates and can accommodate any type of model. Durable constructed and well-engineered, the P-6/3 offers many types of available mounting options, including various table, bench, and chair mounts.



Fig 1. ANA-4



Fig 2. AG-3



Fig 3. A-PZ



Fig 4. P-6/3

Power in Control



Superior cutting efficiency – Safe preparation

- Reliable mechanical glide path
- Efficient dentin removal and precise shaping
- No unnecessary loss of tooth substance

All you need to know in your free M_{two}[®] User Card:
www.vdw-dental.com

m_{two}[®]

The efficient NiTi system



Endo Easy Efficient[®]

phone: 0800 30 10 10

fax: 0800 486 9336

www.gunz.co.nz