

PROSTHODONTIC

Replacement of missing teeth and related mouth or jaw structures with bridges, dentures or other artificial devices is an increasingly important part of daily surgery activity.

The impression-taking procedure is the key to success and we supply a complete range of materials to help you create the perfect impression everytime. From bite registration and impression materials through to cements and tissue management, we have got you covered.

For all your Prosthodontic needs Talk To Us.



SECTION 13

PROSTHODONTIC

CONTENTS

MIXING BOWLS	227
ALGINATES	227
CEMENTS - TEMPORARY	227
BITE REGISTRATION MATERIAL	227
TEMPORARY FILLING MATERIAL	227
CEMENTS - PERMANENT	227
TEMPORARY CROWN & BRIDGE	228
IMPRESSION MATERIAL	228-229
ANTI BACTERIAL SCRUBS	230
OXYGEN INHIBITOR	230
IMPRESSION TRAY ADHESIVES	230
IMPRESSION TRAYS	231-232
IMPRESSION TRAY CLEANER	232
RELINING MATERIAL	233
HEMOSTATIC TISSUE MANAGEMENT	233-234
RETRACTION CORD	235
CROWN FORMS	
FRASACO	236-237
PD	238-239
TRANSPARENT STRIP CROWNS	239
PEDODONTAL STRIP CROWNS	240
MIXING DEVICES, CANNULAS & TIPS	240-241

Ask your Gunz Dental representative about the Kettenbach Mixing Machine.

MIXING BOWLS



Bowls
AINSWORTH

Self cleaning silicone mixing bowls dark green in colour.

Flex-O-Bowl	Medium
AIN304321	
Flex-O-Bowl	Large
AIN304322	
Flex-O-Bowl	Extra Large
AIN304323	

ALGINATES



Kromopan
LASCOD

A stable Class A, dust free, thixotropic, fast set alginate with chromatic phase indicator and physical working time adjustment. The stability of the hydrocolloid within Kromopan makes it possible to pour the impression up to 100 hours after impression taking. It is highly compatible with plaster, reproducing all details.

Bulk Pack 20 x 450 gram carton
LASKRM302
450 gram pack
LASKRM302Z

CEMENTS TEMPORARY



Provicol QM
VOCO

Eugenol-free temporary luting cement with calcium hydroxide. For temporary luting of inlays, onlays, crowns, partial crowns, bridges, attachments and posts or temporary filling of small cavities.

QuickMix Syringe 10m., Mixing tips type 10
VOC1074

BITE REGISTRATION MATERIAL



Futar
KETTENBACH

Futar – the convenient classic: hard A-Silicon with convenient working time, accurate results and the utmost comfort.

Futar Bite Registration, 2 x 50ml cartridge + 12 mixing tips

KET11912



Futar D
KETTENBACH

Futar D – extra hard bite registration material with convenient working time, accurate results, no bounce.

Futar D Bite Registration, 2 x 50ml cartridge + 12 mixing tips

KET11932

Futar D FAST Bite Registration, 2 x 50ml cartridge + 12 mixing tips

KET11961



Registrado X-tra
VOCO

A highly accurate bite registration material based on A-silicone technology. Featuring highly

thixotropic and grinding properties, Registrado X-tra provides very precise impressions of occlusal bite registration, fixing facebow registration and retaining gothic arch registrations. It also offers optimised final hardness and a very short setting time of 40s.

Registrado X-tra 2x50ml safety cartridge
VOC2035
Mixing Tips Type 8
VOC2196

TEMPORARY FILLING MATERIAL



Clip F
VOCO

A fluoride releasing Clip. Suitable for temporary fillings, temporary sealing of cavities, especially in inlay/onlay technique.

Syringes 2x 4 g
VOC1289
Tripack Syringes 2x 4 g
VOC1283

CEMENTS PERMANENT



Implant Cement
PREMIER

Ideal for long-term temporaries and any implant retained crowns that may require adjustment or re-treatment. With no chemical bonding, crowns can be removed more easily – if needed. Premier Implant Cement will not wash out!

Contains: 5ml Automix Dual Barrel Syringe, 10 Mixing Tips, 1 x Mixing Pad
PRE1450
Contains: 25 Automix Dual Barrel Syringe Mixing Tips
PRE1451

TEMPORARY CROWN & BRIDGE



Structur Premium

VOCO

Universal use for highly aesthetic temporary crowns and bridges, inlays and onlays for the most ambitious demands.

Intro Kit Contains: 75gm cartridge A3, dispenser type 2, & mixing tips type 6
VOC1710

Refills

76 gm cartridges & 8 mixing tips

- A1 VOC1712
- A2 VOC1713
- A3 VOC1714

Mixing Tips type 6 – Pack of 50

VOC2185



Structur 2 SC

VOCO

A self-curing temporary crown and bridge material, and with properties such as low shrinkage results in an excellent fit and a tension-free material. High material hardness and transverse strength allow a strong functional loading of the temporary piece and due to its minimal water uptake, Structur 2 will stay dimensionally stable in the moist oral environment and does not discolour. Therefore, it is well suited for a long term temporary.

Intro Kit Contains: 75gm cartridge A3, dispenser type 2, & mixing tips type 6

VOC1478

Refills

75 gm cartridges & 8 mixing tips

- A1 VOC1479
- A2 VOC1480
- A3 VOC1481

Mixing Tips type 6 – Pack of 50

VOC2185

IMPRESSION MATERIAL



Monopren Transfer

KETTENBACH

Monopren® Transfer, the A-silicone for uncompromisingly accurate impressions, is a medium viscosity impression material for the monophasic impression technique. With its initial hydrophilicity, high final hardness and optimal elasticity, it is particularly suited to fixation impression for combined prosthetic work and for implantology.

Monopren transfer: 2 x 50ml cartridge + 12 mixing tips

KET11846

Monopren transfer Intro Pack:

1 x 362ml foil bag, 10 mixing tips, 1 cartridge body

KET29096

Monopren transfer Refill pack: 2 x 362ml foil bag

KET11881

Panasil Binetics Putty

KETTENBACH

Panasil® Binetics putty is the first genuine putty in a foil bag. Pure comfort in a unique form. Intelligent benefits: no more manual kneading, exact dispensing and reduced material use. With the convenient Renfert Sympress mixing machine Panasil® binetics putty can be mixed cleanly, safely and quickly in electrical mixing and dosing machines.

Panasil Binetics Putty Soft Intro Kit :

1 x 362ml foil bag, 10 dynamic mixers, 1 cartridge body

KET28104

Panasil Binetics Putty Soft 2x362ml foil bags

KET11671



Panasil Tray

KETTENBACH

Panasil Tray is a high viscosity, additional-curing precision impression tray material, ideal for two-step impression technique (in combination with Panasil initial contact x-light) Panasil Tray Fast heavy has the



IMPRESSION MATERIAL

consistency of a "Heavy Body" with a high final hardness, is easy to cut and offers the basis for efficient work.

Available in 1:1 cartridge (50ml) for easy handling in the manual dispensing gun or in the convenient Renfert Sympress System in a 5:1 foil bag (362ml)

Panasil Tray Fast Heavy Intro Kit: 1 x 362ml foil bag, 10 dynamic mixers, 1 cartridge body

KET28098 

Panasil Tray Fast Heavy 2x362ml foil bags

KET11621 

Panasil Tray Soft Heavy Intro Kit 1x362ml,

KET28102 

Panasil Tray Soft Heavy 2x362ml foil bags

KET11641 



Panasil Initial Contact

KETTENBACH

The low viscosity or extra low viscosity impression materials in the Panasil® family are ideally coordinated to the Panasil® tray materials and offer precise results in the two-step impression technique or for application in the double mix technique.

Panasil initial contact 2x50ml

KET13411 

Ask your Gunz Dental representative about the Kettenbach Mixing Machine.



**Kellys
PREVIDENT**

60 gm base, 60 gm hardener & mixing pad

PVDWP11350 



Identium

KETTENBACH

Identium® is a Vinylsiloxanether® (VSXE®) which combines for the first time the advantages of A-silicone and polyether in one. A perfect balance of the new Vinylsiloxanether® (VSXE®) means the material is extremely hydrophilic, highly flowable, provides excellent resilience and offers a high degree of comfort for the user and the patient. No bad tastes or odours.

Identium® is available in the viscosities Heavy, Medium and Light.

Identium Heavy Intro Kit

1 x 362ml foil bag, 10 dynamic mixers, 1 cartridge body

KET28237 

Identium Heavy Refill Pack

KET13776 

Identium Medium Intro Kit

1 x 362ml foil bag, 10 dynamic mixers, 1 cartridge body

KET28233 

Identium Medium Refill Kit

KET13736 



Identium Light

KETTENBACH

Identium® Light is a low viscosity impression material available in the 1:1 cartridge. It is particularly suitable for the double mix technique (in combination with high viscosity Identium® Heavy).

Identium Light Refill 2x50ml

KET13701 

ANTIBACTERIAL SCRUBS



Consepsis
ULTRADENT

Use prior to cementation, luting (provisional and/or permanent) and direct restorative placement. Consepsis helps address micro-organisms in dentinal tubules. Minimise post op pulpitis and sensitivity by thoroughly cleaning and disinfecting the preparation prior to sealing and restoring.

Consepsis Kit

Contains: 4x 1.2ml Syringe, 20 Blue Mini Dento Infusor Tips & 20 Black Mini Brush Tips

ULT0490

Consepsis Refill - Contains: 20 x 1.2ml Syringes

ULT0491

OXYGEN INHIBITOR



DeOx®
ULTRADENT

A viscous, glycerine-based gel designed to prevent formation of an oxygen inhibition layer on the surface of resin materials when they are polymerised.

Kit 4 x 1.2ml syringes 20 Black mini tips

ULT0238

IMPRESSION TRAY ADHESIVES



Panasil Tray Adhesive
KETTENBACH

For the use with all PVS materials.

Panasil Adhesive 10ml

KET14101



Identium Tray Adhesive
KETTENBACH

For the use with Identium Tray Material (Vinylsiloxanether).

Identium Adhesive 10ml

KET14204



Attach 2 Alginate
ADM

Designed to dry quickly when sprayed on impression trays, so that cohesive failure of Attach2 is minimised and the alginate can be placed in the tray almost immediately. No toxic aromatic compounds are used as solvents

Alginate tray adhesive in a 215 mL aerosol can.

ADM-ATT2C

Alginate tray adhesive in a 60ml bottle with brush within cap.

ADM-ATT2B



VectorVP & BlackMagic
ADM

Vector VP - Vinyl polysiloxane tray adhesive. Silicone tray adhesive in 60 ml bottle, brush in cap
ADM-VP60

BlackMagic - Polyether tray adhesive. Polyether tray adhesive in a 60 mL bottle with brush in cap.

ADM-BMG60



DZolve
ADM

DZolve can be used to thin down tray adhesives that have thickened.

Universal tray adhesive thinner in a 100 mL bottle

DZLV100



Hold
PREVIDENT

For alginate and hydrocolloid impressions.

80ml spray can

PVDWP11461

14ml liquid

PVDWP11460



Brushtix Disposable Tray Adhesive Brushes
ADM

Disposable single use brushes to be used with all ADM tray adhesives. The brushes can be removed from the caps by simply twisting them off.

Single use adhesive brushes. 20 in each packet. Suitable for VP60, ATT260 and BMG60 adhesives.

ADM-BS20

IMPRESSION TRAYS



Triple Tray®
PREMIER



Proven accurate in over 80 million impressions, these trays have a very thin, loose webbing design that will not distort when the patient closes. These trays are suitable for the dual-arch technique.



Posterior
PRE6203

Box 48



Anterior
PRE6213

Box 35



Quadrant
PRE6223

Box 40



Sideless
PRE6243

Box 48



Extended Quadrant
PRE6253

Box 35



3 Quarter Arch
PRE6273

Box 35



Wide body Posterior

19% longer and 33% wider than the original posterior tray. It can capture details of up to 6 teeth.

Box 48

PRE6283



Track It Trays

PROGRESS

These disposable single arch impression trays have super smooth edges to avoid tissue damage, and the lower wall design ensures patient comfort. The posterior palatal lip helps impede impression material overflow and patient gagging. Available in packs of 12.



#1 Full Arch Upper Large

PRO0112



#2 Full Arch Lower Large

PRO0212



#3 Full Arch Upper Medium

PRO0312



#4 Full Arch Lower Medium

PRO0412



#5 Full Arch Upper Small

PRO0512



#6 Full Arch Lower Small

PRO0612



#7 Quadrant Upper Left/Lower Right

PRO0712



#8 Quadrant Upper Right/Lower Left

PRO0812



#9 Upper Anterior

PRO0912



Light Cured Tray Material

W&P



Plaque photo light cured tray blanks.

Pink - Pack of 50

WP-5001



Blue - Pack 50

WP-5016



Natural - Pack 50

WP-5034



IMPRESSION TRAYS



**Master
PREVENT**

Available in packs of 12.

No. 1 Large Upper
PVD11211

No. 2 Large Lower
PVD11212

No. 3 Medium Upper
PVD11213

No. 4 Medium Lower
PVD11214

No. 9 Anterior
PVD11219

No. 5 Small Upper
PVD11215

No. 6 Small Lower
PVD11216

No. 7 Partial Upper Left/Lower Right
PVD11217

No. 8 Partial Upper Right/Lower Left
PVD11218



SWORD

SUP B1 Upper
SW04972/111

SUO B3 Upper
SW04972/113

SS INF B1 Lower
SW04972/211



**Exploit Trays
LASCOD**

Set of 6 perforated satin finish, one piece impression trays. Each box contains small, medium and large.

LASEXK603

Lower PerforMatFinish S/S Large x1

LASEXF215

Lower PerforMatFinish S/S Medium x1

LASEXF214

Lower PerforMatFinish S/S Small x1

LASEXF213

Lower PerforMatFinish S/S Large x1

LASEXF115

Lower PerforMatFinish S/S Medium x1

LASEXF114

Lower PerforMatFinish S/S Small x1

LASEXF113

IMPRESSION TRAY CLEANER



**All Cleaner
LASCOD**

A cleaning fluid for the removal of alginates, plaster and cements from stainless steel, plastic and other alloys instruments and trays.

Lascod All Cleaner Concentrate 1 Litre

LASALC100



**Traypurool
VOCO**

A proven biodegradable cleaning concentrate for the removal of alginates, zinc oxides, cements and adhesive varnishes from impression trays and instruments.

1 litre
VOC2289

IMPRESSION TRAY CLEANER



**NoAL Alginate Tray Cleaner
ADM**

Actively removes all alginates, tray adhesives, zinc oxides and most cements from impression trays and instruments. It is non foaming and perfect for ultrasonic cleaners.

1 litre alginate tray cleaner - 1:10 dilution

ADM-NoAL1

5 litre alginate tray cleaner - 1:10 dilution

ADM-NoAL5

Dental tray adhesive container with integrated basket. 3.5 litre capacity

ADM-VAT



**Detach Dental Tray Adhesive Remover
ADM**

DETACH will remove polyether, polysulphide, polyvinyl siloxane and alginate tray adhesives. DETACH comes in a handy 100ml bottle and a 215ml Aerosol and is used by simply applying a small amount to a piece of gauze and gently rubbing the impression tray.

Dental tray adhesive remover in a 100mL bottle.
ADM-DTCH 100

ADM-DTCH 100

Dental tray adhesive remover in a 215mL aerosol
ADM-DTCH 215

ADM-DTCH 215



**Solvitan
VOCO**

Cleaning agent with citrus fruit aroma, suitable for the removal of paste, silicone or polyether-based impression material from hands and instruments. Suitable for surgeries and labs.

250ml bottle

VOC2312

RELINIE MATERIAL



Ufi Gel Hard C
VOCO

A direct hard relining material perfectly mixed from the cartridge for simple, time saving and bubble-free application. It is methylmethacrylate-free and has a neutral taste and colour. It is suitable for direct and indirect, total or partial relining, lengthening of denture rims and the repair and enlargement of dentures.

Intro Set Contains: 80g cartridge, dispenser type 2, 10 ml adhesive & mixing tips type 8

VOC2218

Set Contains: 80gm cartridge, 10ml adhesive, & mixing tips type 8

VOC2215

Cartridge Refill Contains: 80gm cartridge, & mixing tips type 8

VOC2216

Refill Adhesive Contains: 10ml adhesive

VOC2217



Ufi Gel SC
VOCO

A permanent soft A silicone relining material for full or partial dentures and implants. It can also offers superior cushioning of sharp alveolar processes maximising patient comfort.

Set Contains: 50ml cartridge, 10ml adhesive, 10ml glazing base, 10ml glazing catalyst, mixing tips type 8

VOC2041

Mixing Tips Type 8 – Pack of 50.

VOC2187

HEMOSTATIC TISSUE MANAGEMENT



ViscoStat®
ULTRADENT

ViscoStat Original and Wintermint are viscous 20% Ferric Sulphate coagulative haemostatic gel. ViscoStat initiates profound haemostasis in seconds while demonstrating unsurpassed kindness to tissue and providing elimination of sulcular fluid. ViscoStat® saves chair time and decreases costly impression remakes. Aside from being an effective haemostatic agent, ViscoStat is an excellent alternative medicament to formocresol when performing vital pulpotomies.

ViscoStat Indispense Refill Syringe

ULT0645

ViscoStat Wintermint 30ml Indispense Refill

ULT0643

Wintermint Tissue Management Kit Contains: 30ml IndiSpense syringe, #00, #0, #1 and #2 Ultrapak cord, Ultrapak Packer, Ultrapak cord organiser, 20x Metal Dento-Infusor tips, 10 1.2ml empty Syringes

ULT6401

ViscoStat Econo Kit. 1.2ml Syringes x 20

ULT3070

ViscoStat Mini Kit 1.2ml Syringes x 2

ULT1087



ViscoStat Clear
ULTRADENT

ViscoStat Clear is a new generation viscous and spreadable 25% Aluminium Chloride haemostatic gel. The product is ideal to use in the aesthetic zone and quickly stops minor bleeding and sulcular fluid. The gel leaves no residue and rinses off with ease. The direct syringe delivery system provides added convenience in placement and eliminates waste. It is available in four convenient kit configurations.

Indispense Kit: 1 x 30ml Indispense, 20x Metal Dento-Infusor Tips, 20 x 1.2 ml Syringes

ULT6407

Indispense Refill: 1 x 30 ml Syringe Refill

ULT6408

Syringe Mini Kit: 4 x 1.2 ml prefilled syringes, 20 x Metal Dento-Infusor Tips

ULT6409

Econo Refill: 20 x 1.2 ml Prefilled Syringes

ULT6410



Astringent®
ULTRADENT

The only American Dental Association accepted haemostatic solution available. Astringent® is an aqueous, 15.5% ferric sulphate haemostatic solution that is more tissue kind than other ferric sub-sulphate solutions. It is highly recommended as being effective and easy to use for control of bleeding, sulcular fluid and an excellent alternative to formocresol when performing vital pulpotomies.

30ml bottle

ULT0111

30ml Indispense refill syringe

ULT0686



Astringent®X
ULTRADENT

Astringent X is an aqueous 12.7% iron solution containing equivalent ferric Sulphate and ferric sub-Sulphate for difficult-to-stop, problem bleeding.

30ml Bottle

ULT0112

USE DENTO INFUSOR TIPS WITH ANY ULTRADENT HEMOSTATIC SOLUTION TO ACHIEVE OPTIMUM HEMOSTASIS. SEE IMPLANT SECTION PAGE 131.

HEMOSTATIC TISSUE MANAGEMENT



Hemodent™
PREMIER

Effectively stops minor gingival bleeding. Epinephrine free to help avoid cardiac reactions. Recommended to arrest bleeding during crown preparations, before impression taking, placing of fillings, cement bases or cementing inlays or crowns.

10cc Bottle

PRE7071

40cc Bottle

PRE7073



Hemostasyl
PIERRE ROLLAND

Haemostatic dressing with mechanical action to stop moderate bleeding and sulcular fluid in all routine dental procedures

Thanks to the high thixotropicity and astringent properties of 15% Aluminium Chloride, Hemostasyl achieves efficient haemostasis in less than 2 minutes by combining both chemical and mechanical actions when applied with a calculate and moderate pressure.

The disposable angled cannulas provide precise application directly to the bleeding reducing the risk of waste and shortening the contact time. Hemostasyl ensures excellent adhesion to the gingival mucosa without any compression and can be washed away with a spray afterwards. The paste is completely painless on healthy periodontium.

Starter Kit

Contains 2x2g syringes, 40 disposable tips

PR-HEMOSKIT

Refills - Pack 2x2g syringes

PR-HEMOSSYR

Pack 40 disposable tips

PR-HEMOSTIPS



Place Traxodent™



After 2 minutes,
rinse away

Traxodent™
PREMIER

Traxodent – Predictable, Absorbent, Fast! Superior haemostasis and gentle retraction in just 2 minutes.

- Fast & Predictable haemostatic paste with 15% aluminium chloride
- Sleek syringes and tips make application fast & easy!
- Replaces 2nd cord in double cord technique
- Economical with low ongoing costs



Premier Traxodent Value Pack x25 syringes and 50 applicator tips

PRE7091

Premier Traxodent Starter Pack x7 syringes and 15 applicator tips

PRE7093

Premier Traxodent Applicator Tips X50

PRE7092



Expasyl™

PIERRE ROLLAND/ACTEON GROUP

Expasyl is injected into the sulcus by using an applicator. A sulcular opening is obtained in between 30 seconds and 2 minutes, depending on gum texture. The end result is a dry, wide open sulcus with access to the cervical margins.

Expasyl is suitable for use with dental impressions with conservative dentistry, when restoring class II and V caries, placement of rubber dam, placement of orthodontic bands, and apical surgery.

Expasyl Premium Kit

Contains 1 x Applicator, 10 x Expasyl carpules (strawberry), 40 x single use tips

PR-EXPKPRE

Refills

Pack 20 refill strawberry carpules

PR-EXCAP20

Pack 10 refill strawberry carpules

PR-EXCAP10

Pack 20 refill original carpules

PR-EXPR

Pack 40 disposable tips

PR-EXPC

Pack 100 disposable tips

PR-EXPTIP100

Pack 40 disposable tips, pre bent

PR-EXPTIPBENT40

Pack 100 disposable tips, pre bent

PR-EXPTIPBENT100

Expasyl Dispenser

Can be autoclaved

PR-260900



RETRACTION CORD



Ultrapak® Retraction Cord
ULTRADENT

Ultrapak cord is made of 100% cotton, knitted into thousands of tiny loops forming long interlocking chains. The unique knitted design exerts a gentle continuous outward force following placement as knitted loops seek to open. Ultrapak cord is designed to enhance the ViscoStat/Astringedent and Dento-Infusor Tissue Management Technique; however, conventional techniques using alum, aluminum chloride, etc., are also enhanced when using Ultrapak plain knitted cords, as they carry significantly greater quantities of haemostatic solution than conventional cords.

Ultrapak cord size 0

ULT0131 

Ultrapak cord size 1

ULT0132 

Ultrapak cord size 2

ULT0133 

Ultrapak cord size 3

ULT0134 

Ultrapak cord size 00

ULT0136 

Ultrapak cord size 000

ULT0137 

UltraPak Kit (Sizes: 00, 0, 1, 2.)

ULT0130 



Fisher's Ultrapak® Packers
ULTRADENT

These specially designed packers enhance the ease and completeness of packing Ultrapak knitted cord. The slim, slightly serrated heads allow the packer to easily slip the cord into the gingival sulcus.

Ultrapak Packer Regular

ULT0170 

Ultrapak Packer Small

ULT0171 



Hemodent Cord
PREMIER

Hemodent Cords already are impregnated with Hemodent Liquid. The cord has a pliable consistency to be easily packed subgingivally without fraying.

Braided Cord Medium – Thin

PRE7067 

Twisted Cord #3

PRE7077 

Twisted Cord #9

PRE7078 

No matter
what your
dental needs
we've
got you
covered...

- ★ Clinical
- ★ Equipment
- ★ Laboratory
- ★ Implant
- ★ Orthodontics

TALK TO US

FRASACO - CROWN FORMS

FOR PEDODONTAL SIZES REFER TO PAGE 240.

STRIP CROWN

FRASACO



Can be used for temporary restoration of crown preparations and as a matrix for fractured corners and end restorations. 5 crowns per refill.

CENTRAL INCISORS - RIGHT UPPER

116	115	114	113	112	111	110
#116	FRA116-5					
#115	FRA115-5					
#114	FRA114-5					
#113	FRA113-5					
#112	FRA112-5					
#111	FRA111-5					
#110	FRA110-5					

CENTRAL INCISORS - LEFT UPPER

210	211	212	213	214	215	216
#210	FRA210-5					
#211	FRA211-5					
#212	FRA212-5					
#213	FRA213-5					
#214	FRA214-5					
#215	FRA215-5					
#216	FRA216-5					

LATERAL INCISORS - RIGHT UPPER

125	124	123	122	121	120
#125	FRA125-5				
#124	FRA124-5				
#123	FRA123-5				
#122	FRA122-5				
#121	FRA121-5				
#120	FRA120-5				

LATERAL INCISORS - LEFT UPPER

220	221	222	223	224	225
#220	FRA220-5				
#221	FRA221-5				
#222	FRA222-5				
#223	FRA223-5				
#224	FRA224-5				
#225	FRA225-5				

CANINES - RIGHT UPPER

136	135	134	133	132	131	130
#136	FRA136-5					
#135	FRA135-5					
#134	FRA134-5					
#133	FRA133-5					
#132	FRA132-5					
#131	FRA131-5					
#130	FRA130-5					

CANINES - LEFT UPPER

230	231	232	233	234	235	236
#230	FRA230-5					
#231	FRA231-5					
#232	FRA232-5					
#233	FRA233-5					
#234	FRA234-5					
#235	FRA235-5					
#236	FRA236-5					

FIRST PREMOLARS - RIGHT UPPER

144	143	142	141
#144	FRA144-5		
#143	FRA143-5		
#142	FRA142-5		
#141	FRA141-5		

FIRST PREMOLARS - LEFT UPPER

241	242	243	244
#241	FRA241-5		
#242	FRA242-5		
#243	FRA243-5		
#244	FRA244-5		

FRASACO - CROWN FORMS

FOR PEDODONTAL SIZES REFER TO PAGE 240.

FIRST MOLARS - RIGHT UPPER



163	162	161	160
#163	FRA163-5		
#162	FRA162-5		
#161	FRA161-5		
#160	FRA160-5		

FIRST MOLARS - LEFT UPPER



260	261	262	263
#260	FRA260-5		
#261	FRA261-5		
#262	FRA262-5		
#263	FRA263-5		

FIRST MOLARS - RIGHT LOWER



463	462	461	460
#463	FRA463-5		
#462	FRA462-5		
#461	FRA461-5		
#460	FRA460-5		

FIRST MOLARS - LEFT LOWER



360	361	362	363
#360	FRA360-5		
#361	FRA361-5		
#362	FRA362-5		
#363	FRA362-5		

FIRST PREMOLARS - RIGHT LOWER



444	443	442	441
#444	FRA444-5		
#443	FRA443-5		
#442	FRA442-5		
#441	FRA441-5		

FIRST PREMOLARS - LEFT LOWER



341	342	343	344
#341	FRA341-5		
#342	FRA342-5		
#343	FRA343-5		
#344	FRA344-5		

CANINES - RIGHT LOWER



434	433	432	431	430
#434	FRA434-5			
#433	FRA433-5			
#432	FRA432-5			
#431	FRA431-5			
#430	FRA430-5			

CANINES - LEFT LOWER



330	331	332	333	334
#330	FRA330-5			
#331	FRA331-5			
#332	FRA332-5			
#333	FRA333-5			
#334	FRA334-5			

LATERAL INCISORS

RIGHT LOWER

LEFT LOWER



422	421	420	320	321	322
#422	FRA422-5				
#421	FRA421-5				
#420	FRA420-5				
#320	FRA320-5				
#321	FRA321-5				
#322	FRA322-5				

CENTRAL INCISORS

RIGHT LOWER

LEFT LOWER



413	412	411	311	312	313
#413	FRA413-5				
#412	FRA412-5				
#411	FRA411-5				
#311	FRA311-5				
#312	FRA312-5				
#313	FRA313-5				

PD - PELLA CROWN FORMS

PD - Pella

Made in cellulose acetate for the peel off technique, they reproduce the original form of the tooth. Right crowns are identified in yellow and left crowns in pink. All packs contain pairs – left and right.

UPPER CENTRALS – NORMAL

R1+L1	PDCF1	
R11+L11	PDCF11	
R21+L21	PDCF21	
R31+L31	PDCF31	
R41+L41	PDCF41	
R51+L51	PDCF51	

EXTRA LARGE CENTRALS

R161+L161	PDCF161	
R171+L171	PDCF171	

UPPER CANINES

R13+L13	PDCF13	
R3+L3	PDCF3	
R23+L23	PDCF23	

UPPER CENTRALS – CERVICAL

R101+L101	PDCF101	
R111+L111	PDCF111	
R121+L121	PDCF121	
R131+L131	PDCF131	
R141+L141	PDCF141	
R151+L151	PDCF151	

LOWER CANINES

R203+L203	PDCF203	
R213+L213	PDCF213	

EXTRA LARGE CANINES

R33+L33	PDCF33	
R43+L43	PDCF43	

UPPER LATERALS

R12+L12	PDCF12	
R2+L2	PDCF2	
R22+L22	PDCF22	
R32+L32	PDCF32	
R42+L42	PDCF42	

LOWER INCISORS

R211+L211	PDCF211	
R212+L212	PDCF212	
R201+L201	PDCF201	
R202+L202	PDCF202	


PD - CROWN FORMS

BICUSPIDS



R14+L14	PDCF14	
R15+L15	PDCF15	
R24+L24	PDCF24	
R4+L4	PDCF4	
R5+L5	PDCF5	

MOLARS



Upper R16+L 16	PDCF16	
Upper R17+L 17	PDCF17	
Lower R6+L 6	PDCF6	
Lower R7+L 7	PDCF7	

ASSORTED

Assortment L&R 39	PDCFASS	
-------------------	---------	--

TRANSPARENT STRIP CROWNS



Anterior assortment



Posterior assortment

Anterior Assortment

FRASACO

174 transparent strip crowns in plastic box with transparent cover. Three each of every size in 58 compartments.

FRAST-FK174	
-------------	--

Posterior Assortment

FRASACO

96 transparent strip crowns in plastic box with transparent cover. Three each of every size in 32 compartments.

FRAST-BK96	
------------	--



Compact Assortment

FRASACO

Current assortment (current anteriors and posteriors). 58 transparent strip crowns in plastic box with transparent cover. One strip crown of each size in 58 compartments.

FRAST-FBK58	
-------------	--

PEDODONTAL STRIP CROWNS

FRASACO

Placing prefabricated Frasaco pedodontal strip crowns is an excellent technique for retaining badly damaged deciduous teeth. Previously, deciduous teeth with carious defects were extracted once these exceeded a certain stage. Today they are left in the mouth as a spacer as it is now known that the deciduous teeth are important for the unimpeded development of children's jaws. The advances made in conservative dentistry have changed the manner in which this field is approached.

CENTRAL INCISORS

514	513	512	511	611	612	613	614
#514	FRA514-5						
#513	FRA513-5						
#512	FRA512-5						
#511	FRA511-5						
#611	FRA611-5						
#612	FRA612-5						
#613	FRA613-5						
#614	FRA614-5						

LATERAL INCISORS

524	523	522	521	621	622	623	624
#524	FRA524-5						
#523	FRA523-5						
#522	FRA522-5						
#521	FRA521-5						
#621	FRA621-5						
#622	FRA622-5						
#623	FRA623-5						
#624	FRA624-5						

CANINES

534	533	532	531	631	632	633	634
#534	FRA534-5						
#533	FRA533-5						
#532	FRA532-5						
#531	FRA531-5						
#631	FRA631-5						
#632	FRA632-5						
#633	FRA633-5						
#634	FRA634-5						

MIXING DEVICES & TIPS

Mixing Devices & Tips

KETTENBACH



Applyfix Dispensing Gun (1:1 cartridges)
KET17203



Intraoral Tips White x96
KET17221

Intraoral Tips Yellow x 96
KET17222



Mixing Tips Green x100
KET17235

Mixing Tips Yellow x100
KET17240



Dynamic Mixer Tips x50
KET17253

Ask your Gunz Dental representative about the Kettenbach Mixing Machine.

MIXING DEVICES, CANNULAS & TIPS REFERENCE

Dispenser – Type 2



Dispenser to extrude dental material from cartridges. For 1:1 or 1:2 cartridges (50ml)

Suitable for

Registrado X-tra, Ufi Gel SC, Ufi Gel hard C, Structur 2 SC, Structur 2 SC Long Setting, Structur Premium, Provicol C, Contrast, Voco Pac C.
(1x Dispenser Type 2, complete)

VOC2158



Dispenser – Type 3



Dispenser to extrude dental material from cartridges. For 1:1 or 1:2 cartridges (25ml)

Suitable for

Rebilda SC, Rebilda DC.
(1x Dispenser Type 3, complete)

VOC2259



Applicator for AC capsules



Forceps-type applicator to directly apply the content of Voco application capsules (AC) into the cavity.

(1x AC Applicator)

VOC2331



Activator for AC capsules



Device to activate Voco application capsules (AC).

(1x AC Activator)

VOC9300



Application Cannulas



Type 40 for Fissurit F/FX. Pack of 100.

VOC2145



Type 42 for Vococid. Pack of 30.

VOC2310



Type 44 for Admira Flow, Arabesk Flow, Block Out Gel LC. Pack of 100.

VOC2144



Type 41 for Ionoseal, Calcimol LC, Amaris® Flow, Grandio® Flow, Amaris® Gingiva, Bifix QM Try-in. Pack of 100.

VOC2146



Type 43 for Fluoridin N 5 cartridges. Pack of 100.

VOC2167



Type 45 for Grandio® Seal, Cimara Opaker LC. Pack of 100.

VOC2147



Mixing Tips for Syringes



Type 6 for Provicol C, Structur 2 SC, Structur 2 SC Long Setting, Structur Premium. Pack of 50.

VOC2185



Type 8 for Ufi Gel SC, Ufi Gel Hard C, Voco Pac C. Pack of 50.

VOC2187



Type 10 for Provicol QM, Structur 2 QM, Structur Premium QM. Pack of 50.

VOC2189



Type 12 for Rebilda DC in cartridge, Rebilda SC. Pack of 50.

VOC2192



Type 7 for Contrast medium, Contrast light. Pack of 50.

VOC2186



Type 9 for Fit Test C & B. Pack of 50.

VOC2188



Type 11 for Rebilda DC in the QuickMix syringe, Bifix QM. Pack of 50.

VOC2191



Type 14 for Bifix SE. Pack of 50.

VOC2194



Type 16 (for intraoral tips type 2). Suitable for Registrado X-tra. Pack of 50.

VOC2196



VALO™

broadband LED curing light



“Light Done Right”

- 3 power levels – Standard/High/Plasma to cure all types of material
- Focused collimated Beam for complete, uniform cure
- Lightweight uni-body construction provides unprecedented access and durability

Contact your Gunz Dental representative to arrange your in-surgery demonstration.

 **ULTRADENT**
PRODUCTS, INC.

phone: 0800 30 10 10

fax: 0800 486 9336

www.gunz.co.nz