

# RESTORATIVE

As the most important component of a surgery's daily activities, having a full array of restorative materials on hand is essential. Our range includes everything you need to ensure a successful procedure. Whether you're in the market for an amalgamator, composite, etchant or reinforcement material, we have the most appropriate solutions for you.

For all your Restorative needs Talk To Us.



## SECTION 14

# RESTORATIVE

## CONTENTS

AMALGAMATORS	245
CURING LIGHTS	246
AMALGAM GUNS	247
ACRYLIC SPRUES	247
PINS & POSTS	247-248
REINFORCEMENT MATERIALS	249
BRUSHES	250
CARIES INDICATOR	250
ETCHANT	250
BONDING & ADHESIVE	251
BASE LINER	251-252
CEMENTATION & LUTING	252
FLOWABLE COMPOSITE	252-253
UNIVERSAL COMPOSITE	253-256
POSTERIOR COMPOSITE	256
TEMPORARY FILLING	257
GLASS IONOMER	257
CORE BUILD UP MATERIAL	258
MATRICES & STRIPS	259-260
WEDGES	260
ARTICULATING PRODUCTS	260-261
FINISHING & POLISHING	261-262
COMPOSITE WETTING AGENTS	262
COMPOSITE SEALER & VARNISH	263
DESENSITISER	263

## AMALGAMATORS

Softly  
SATELEC

A universal high speed amalgamator for pre-measured amalgam capsules. This 230V unit is simple and easy to use. It is microprocessor controlled to ensure repeatability of selected time. The oscillating range of 25mm and a mixing frequency of 4,000 vibrations per minute provides a perfect homogenous mix and an ideal consistency for the type of material that you choose.

To guarantee safety, the unit is deactivated once the front cover is opened to access the capsules. It is stable, extremely quiet and, being made of injection moulded plastic, is easy to clean and maintain.

Amalgamator (210 x 205 X 120mm)

SATF03160



## CURING LIGHTS



**VALO**  
ULTRADENT'S NEW LED CURING LIGHT

## KIT

## VALO Kit

- 1 - VALO LED curing light - 7' cord
- 1 - Power supply with universal plugs - 6' cord
- 1 - Handpiece bracket holder
- 50 - VALO Barrier sleeves

ULT5919



## BARRIER SLEEVE

## Barrier Sleeve Refill 500pk

ULT5932



## Valo LED Curing Wand

ULTRADENT

VALO is the LED curing light for unsurpassed accessibility, accuracy and efficiency, all contained within a sleek body and designed to deliver the right power in the right place.

VALO is anything but your ordinary curing light. Its optimised power output, broad spectrum custom LEDs and precisely collimated beam provide consistent and uniform cures. With three curing options (standard power, high power and a new plasma mode), VALO has the ability to cure all types of restorations with the touch of a button. With a three-second time cycle and built-in two-second safety delay, VALO's plasma mode has been said to rival plasma arc lights in its level of rapid cure.

No longer do you have to decide if you want a normal output, a high intensity output or the power of a PAC light. VALO has all of that in one light, and in a very light and compact unit. Its ability to produce the greatest amount of energy at distances up to 10 mm from the tooth is just amazing. VALO's sleek, aluminum aerospace wand-style body is lightweight and exceptionally durable.

- 3 power levels – Standard/High/Plasma to cure all types of material
- Focused collimated Beam for complete, uniform cure
- Lightweight uni-body construction provides unprecedented access and durability

VALO Curing Wand Kit 230V

ULT5919

VALO Barrier Sleeves Pkt 500

ULT5932

CURING LIGHTS



**Satelec Mini LED Ortho 2**

Supplied with a Ø 5,5mm light guide and comes standard with same charger as the MiniLed AF.

- Features a large LCD screen
- It offers 4, 8, 12 and 32 second curing cycles.
- High power of up to 3,000 W/cm<sup>2</sup> using a Ø 5,5mm light guide
- Wave length is 420-480nm to cover all different composites but also sealants and bonding agents
- It prompts the user by offering four different alert modes – beep, micro flash, micro flash and beep or no additional alert
- Standby mode is active after 5 minutes
- The Lithium-Ion battery can be re-charged at any time and will last to cordlessly bond a minimum of 15 patients
- Exclusive patented heat dissipation method prevents over-heating and protects circuitry - no fan needed.

Mini LED Ortho 2

SATF055220 

**Satelec Mini LED Light Guides**

Perfect for all polymerisation of composite.

No 1 Opalescence multi-fibered 7.5mm Guide

SATF02648 

Perfect for ortho brackets, small posterior, anterior and cervical areas.

No 2 Opalescence multi-fibered 5.5mm Guide

SATF02652 

The 90 degree angle is designed to deliver light into hard to reach areas such as root canals and pulpotomies. Perfect for lingual work such as orthodontics & splinting; ideal for paediatric dentistry.

No 3 Amber multi-fibered 3.5mm Guide

SATF02651 



**MINI LED CORDLESS**

**Mini LED Cordless Curing Light**

SATELEC

Cordless curing light with optimal power, for efficiency and speed.

- Emits light in the most relevant and efficient part of the spectrum, to cure all photo initiated products
- Very fast polymerisation (i.e. 6-12 seconds) to cure 2mm of composite, in fast mode
- Choice of 3 modes (fast, pulse and ramp)
- Latest generation of CMS electronics
- Ergonomic compact design weighing only 190g

Curing Light Grey: 3 modes: Fast, Pulse & Ramping. Includes 7.5mm opalescent light guide, battery and protective light shield.

SATF02621 

Mini LED Curing Light Battery (Grey)

SATF02520 

**Mini LED AF Curing Light**

SATELEC

World's first cordless curing light with "Autofocus" function.

- More efficient, faster curing is achieved by the automatic calculation of precise exposure time measured in accordance with the distance
- Better curing achieved by most powerful (power density of 2,200 to 3,000 mW/cm<sup>2</sup>) and highly efficient light spectrum for deep polymerisation of composites
- Three curing modes (fast, pulse and ramping)
- Laser target ring allows for easier and precise positioning of the light guide
- Incorporated radiometer to enables instant check of the power density, Lithium-Ion battery

Mini LED AF Light Guide - opalescent 7.5mm

SATF02654 

Mini LED AF Light Guide - opalescent 5.5mm

SATF02655 


Mini LED AF Curing Guide

SATF02900 

Mini LED AF2 OEM

SATF02980 


## AMALGAM GUNS



PD

Stainless steel and autoclavable, with interchangeable heads.

- Duo RA
- PDAGDRA
- Duo Regular
- PDAGDST
- Duo Gun Spring Only
- PDAGDSPRING
- Duo Gun Interior Stem & Spring
- PDAGDSSSET
- Duo Gun Head Nozzle Only
- PDAGDHNRA
- Duo Gun Straight Head Nozzle Only
- PDAGDHNSTR



PREMIER

Fully Autoclavable with removeable parts for cleaning.

- Straight
- PRE62001
- RA
- PRE62006
- Nose Tip Straight
- PRE62011
- Nose Tip RA
- PRE62013
- Plunger Straight
- PRE62015
- Plunger RA
- PRE62017

## ACRYLIC SPRUES



PD

- PDAS1 Size 1
- PDAS2 2
- PDAS3 3
- PDAS4 4
- PDAS6 6

## PINS & POSTS



### IntegraPost™ Ultra-Passive Integrating Core Retainer



PREMIER



The IntegraPost features an innovative round Flow-Thru-Head and a parallel passive design to optimise the adhesive buffer zone, conserve tooth structure and minimise stress.

Made from bio-compatible titanium alloy the IntegraPost is stronger than resin fiber yet more elastic than stainless steel.

The midflange and shoulder are designed to dissipate stress away from the canal. The safe, non threaded parallel shank minimises potential root fracture. Anti-rotational knurling locks into the cement



### Refills — Packs 10

PREMIER

- Size 3
- PRE1105
- Size 4
- PRE1110
- Size 4.5
- PRE1115
- Size 5
- PRE1125
- Size 6
- PRE1130



### IntegraPost System Mini Intro Kit

PREMIER

Contains: 6 Posts (2 x Sizes 3, 4 & 4.5), 3 Drills (Sizes 3, 4, & 4.5), 1 x IP Carrier

- PRE1001



### Integra Drills

PREMIER

Ideally sized in relation to post diameter, permitting 2.5 times more cement into the post space. This creates a buffer zone to absorb potentially stressful forces and increases retention.

Packs of 3 Drills.

- Size 3
- PRE1205
- Size 4.
- PRE1210
- Size 4.5
- PRE1215
- Size 5
- PRE1225
- Size 6
- PRE1230



### Precision Post System

PRECISION

These posts are made from the high strength corrosion resistant Titanium alloy. All sizes are compatible with the ParaPost® system from Coltene Whaledent. All posts are parallel sided, fluted for maximum retention and vented to reduce pressure and stress during placement. The sizes 4.5 and 5.5 are a Precision innovation for a more exact fit in-between the most popular post sizes.

- Refill 3 — Pack 10
- PROTVP-3
- Refill 4 — Pack 10
- PROTVP-4
- Refill 4.5 — Pack 10
- PROTVP-4.5
- Refill 5 — Pack 10
- PROTVP-5
- Refill 5.5 — Pack 10
- PROTVP-5.5

PINS & POSTS

Refill 6 – Pack 10

PROTVP-6

Refill 7 – Pack 10

PROTVP-7

Drills are made from highest quality corrosion resistant alloy, and are specially designed to have a non-aggressive cutting edge to form safe, concentric post spaces. The shank has a silicone colour ring to indicate size. The ring retains its colour after repeated sterilisation. All post and drill packages have colour code conversion charts to compare Precision drill and post sizes with ParaPost® sizing.

Refill 3 – Pack of 2

PROD-3

Refill 4 – Pack of 2

PROD-4

Refill 4.5 – Pack of 2

PROD-4.5

Refill 5 – Pack of 2

PROD-5



Crown Posts with Sprues

PD

Size 1 PDCPS1

Size 2 PDCPS2

Size 3 PDCPS3

Size 4 PDCPS4



Rebuild Post  
VOCO

Fibre-reinforced composite post. Endodontic post build-up for support and anchoring of coronal restoration in situations with insufficient tooth substance.

- Dentine-like elasticity behaviour, high transverse strength
- High radiopacity (350%Al)
- Translucency like dentine
- Anatomical shape
- Adhesive luting
- Removable

Intro Set: 5x posts each of (ø 1.2 mm, ø 1.5mm, ø 2.0 mm), 1x drill each of (ø 1.2 mm, ø 1.5 mm, ø 2.0 mm)

VOC1770

REBILDA POST - DENTINE-LIKE E-MODULUS

Refills:

5x posts, Rebuild Post 12 (ø 1.2 mm)  
VOC1772

5x posts, Rebuild Post 15 (ø 1.5 mm)  
VOC1773

5x posts, Rebuild Post 20 (ø 2.0 mm)  
VOC1774

1x drill, Rebuild Post 12 (ø 1.2 mm)  
VOC1777

1x drill, Rebuild Post 15 (ø 1.5 mm)  
VOC1778

1x drill, Rebuild Post 20 (ø 2.0 mm)  
VOC1779



Rebuild Post System  
VOCO

Complete set for post-endodontic treatment for placing 15 posts.

- All materials in the set match each other
- Post insertion and core-build up in one step

**Futurabond DC**

- Safe self-cure
- High adhesion without additional etching

**Rebuild Post**

- E-modulus like dentine
- High radiopacity

**Rebuild DC**

- Suitable for post luting
- Cuts like dentine, good physical properties
- Low setting temperature



POSTS ILLUSTRATIONS ARE NOT ACTUAL SIZE

Complete Set: 5x posts each of (ø 1.2 mm, ø 1.5mm, ø 2.0 mm), 1x drill each of (ø 1.2 mm, ø 1.5 mm, ø 2.0 mm), 1x reamer (ø 0.7mm), 5 ml Ceramic Bond, 15x Futurabond DC SingleDose, 10g QuickMix syringe Rebuild DC dentine, accessories.  
VOC1771

REINFORCEMENT MATERIALS



**BioSplint**

PIERRE ROLLAND / ACTEON GROUP


BioSplint® is a multiple use dental reinforcement ribbon. It is used for stabilisation and immobilisation of loose teeth, immediate replacement of avulsed teeth, reinforcement and construction of both provisional and permanent bridges, or maintenance of interdental spaces. BioSplint ribbon is made of polyethylene terephthalate woven fibers which have undergone various treatments to clean the ribbon perfectly and improve adhesion of the composite. One of these treatments involves the use of supercritical CO2 which guarantees complete elimination of any impurities that may result from the manufacturing process. Biosplint splinting ribbons are very supple and perfectly follow the contours of the teeth and dental arch. BioSplint is a simple technique to use. It is fast, reversible, aesthetic and perfectly painless.

BioSplint Ribbon Refills contain:  
8 x ribbons individually wrapped in aluminum bags (width: approximately 2.5mm, length: approximately 50mm, thickness 0.15mm), 8 x "patient traceability" labels

PR-BSREF 

BioSplint Thermocutter  
Especially suited for a clean cutting of the splint. Works from a battery provided with the device.

PR-BSCUT 

BioSplint Accessories includes:  
8 x single use brushes, 3 x probes which can be autoclaved  
PR-BSACC 



**Ribbond®**  
RIBBOND

Ribbond's unique combination of ultra-high strength fibers, enhanced bondability and patented cross-link lock-stitch leno weave makes the strength and fracture toughness of Ribbond reinforced prostheses unsurpassed by other fiber reinforcements. Ribbond prevents fracture failures in dental composites and acrylics. Ribbond prostheses are strong and maintain their strength with continued use.

**Assorted Starter Kit**

RIBBOND

Original Ribbond is preferred for use with applications in which thinness and higher modulus are primary concerns. These applications include periodontal splints, orthodontic retainers, endodontic post and cores or short span bridges.

Ribbond Starter Kit  
Contains Scissors, Manual, refills 2, 3 and 4mm.

RIBWAST 

**Assorted Refill Kits**

RIBBOND

2mm  
RIBWRE2 

3mm  
RIBWRE3 

4mm  
RIBWRE4 

Assorted  
RIBWREAST 

**Ribbond-THM**  
RIBBOND



Made from thinner fibers with a higher thread count, Ribbond-THM has a higher flexural strength than regular Ribbond and is only 0.18 mm thick. Ribbond-THM offers the same crack-stopping leno-weave and high strength fibers of Original Ribbond with greater manageability, greater ease-of-use, and thinner results.

**Assorted THM Starter Kit**

RIBBOND

Contains everything you need except for the composite & acrylic that you already use in your practice.

RIBASTTHM 

**Refill Kits**

RIBBOND

2mm  
RIBRE2THM 

3mm  
RIBRE3THM 

4mm  
RIBRE4THM 


**Orthodontic THM Starter Kit**

RIBBOND

Especially designed for Orthodontic bonded retainers with a 1mm wide ribbon.

RIBKITOTHM 

Refill Kits  
RIBRE0THM 

Ribbond Universal Scissors  
RIBSCISSORS 

BRUSHES



**ICB® Intracoronar Bristle Brush**  
ULTRADENT

The ICB® brush was designed to clean intra-coronally as well as extracoronally before bonding fillings, inlays, onlays and crowns. It is great for cleaning areas with restricted access. It is autoclavable and designed for use with a slow speed handpiece.



- Pack 30
- ULT1076 
- Pack 100
- ULT1077 

CARIES INDICATOR



**Seek®**  
ULTRADENT



Is the only syringe-delivered caries indicator for quick and easy placement. It contains D&C dyes in a glycol base and stains caries red.

- Kit - 4 x 1.2ml syringes & 20 x Black Mini Tips
- ULT0209 
- Refill - 4 x 1.2ml syringes
- ULT0210 



**Sable™ Seek®**  
ULTRADENT

Contains FD&C dyes in a glycol base and stains carious dentine greenish-black.

- Kit - 4 x 1.2ml syringes & 20 x Black Mini Tips
- ULT0233 
- Refill - 4 x 1.2ml syringes
- ULT0234 

Both are used for rapid easy identification of carious dentine. Clean and precise placement can be achieved with the use of Black Mini Brush Tips.

FOR FULL DETAILS ON MINI BRUSH TIPS SEE PAGE 118.

ETCHANT



**Ultra-Etch®**  
ULTRADENT

Viscous 35% Phosphoric acid gel provides optimum viscosity for precise placement and superior control. Ultra Etch is the only self-limiting etchant on the market with an average depth of etch of 1.9 microns with every 15 second etch. The fumed silica base of Ultra Etch means that gel rinses away cleanly leaving no smear layer to interfere with bonding strengths.

- Kit Contains:
- 4 x .2 ml syringe.
- 10 x Inspiral brush tips
- 10 x Blue blue micro tips

ULT0163 

Refill - 4 x syringes

ULT0164 

Econo Kit - 20 x syringes + 20 micro tips

ULT0167 

35% Econo Kit - 20 x syringes

ULT0168 

35% Indispense Kit

ULT0383 

35% Indispense Refill

ULT0685 



**EtchArrest®**  
ULTRADENT

A sodium bicarbonate and calcium source, acid neutralising medium. Viscous enough not to run

on vertical surfaces, yet can still flow into tiny irregularities. Rinses off with water.

- Kit Contains:
- 4 x 1.2ml syringes
- 20 x Black Micro Tips

ULT0625 



**Porcelain Etch**  
ULTRADENT

Is a specially buffered, viscous hydrofluoric acid designed for intra oral or extraoral porcelain etching. It should be used to etch porcelain or old macro filled composite prior to bonding and repair. Suitable when used with veneers and inlays.

Refill - 2 x syringes

ULT0406 



**Silane**  
ULTRADENT

Ultradent Silane is a single component Silane suitable for use following the application of Porcelain Etch.

Refill - 2 x syringes

ULT0410 



**Gel-Etch**  
PREVIDENT

Gel Etch - 12 gram syringe

PVDPR-7012 

**BONDING & ADHESIVE**



**PQ1® Single Syringe Resin**  
ULTRADENT

A unique, patented one of a kind chemistry. A tissue kind primer and bond resin combination using ethyl alcohol as the solvent carrier. It is a 40% filled, fluoride releasing, radiopaque natural filler. Suitable for 95% of all direct restorative needs and can be used on most materials including dentine, enamel, porcelain and composite.

Intro Kit Contains:

- 2 x 1.2ml PQ1 syringes
- 2 x 1.2ml Ultra Etch 35% syringes
- 20 x Blue micro tips
- 40 x Inspiral brush tips

ULT0615 

Refill - 4 x 1.2ml syringes

ULT0641 



**PermaQuick® Two-Part Bonding System**  
ULTRADENT

A tissue friendly bonding system with an ethyl alcohol carrier and neutral pH primer. It contains natural resins to maximise superior penetration into the dentinal tubules. PermaQuick bonding resin is 45% filled, fluoride releasing, radiopaque and is suitable for direct and indirect bonding procedures. PermaQuick will bond to dentine, enamel, cast metal, porcelain, amalgam and composite.

Primer Refill – 4 x 1.2ml syringes

ULT0633 

Bonding Resin – Refill 4 x 1.2ml syringes

ULT0634 



**Admira Bond**  
VOCO

Light-curing one-component Ormocer®-based dentine/enamel bonding agent. It has one component and one layer with high adhesion and long lasting marginal tightness. The single dose delivery offers a simple, quick and hygienic application with no spillage.

1 x 4ml Bottles

VOC2479 

2 x 4ml Bottles

VOC2432 



**FuturaBond DC**  
VOCO

Dual-curing nano-reinforced self-etch dentine / enamel bond. Chemical-curing with additional DC activator in areas inaccessible with a curing light. DC achieves the same adhesive properties as total-etch preparation without separate etching of the tooth substance. Universal application.

1x Bottle 4 ml Liquid 1, 2x Bottle 4 ml Liquid 2 and accessories

VOC1163 



**FuturaBond DC SingleDose**  
VOCO

Futura bond DC in practical SingleDose. Easy, quick and hygienic – one-step activation blister. No additional devices required. No spilling in any position.

50x SingleDose and accessories

VOC1164 

200x SingleDose and accessories

VOC1165 

**BASE LINER**



**Ultrablend® Plus™**  
ULTRADENT

Ultrablend®Plus™ is a unique, light activated, radiopaque liner material containing calcium hydroxide and calcium hydroxyapatite in a urethane dimethacrylate base. It is suitable for use as a thin cavity liner prior to placing either a composite or amalgam. Proven best for pulp capping.

Kit - Contains:

- 2 x 1.2ml dentin syringes
- 2 x 1.2ml opaque white syringes
- 20 x Black Micro Tips
- 20 x Black Mini Tips

ULT0415 

Dentin Colour Refill - Pack 4 x 1.2ml syringes

ULT0416 

Opaque White Refill - Pack 4 x 1.2ml syringes

ULT0417 



THE DENTAL  
ADVISOR  
+++

**Ionoseal**  
VOCO

Light-cured ready-for-use Glass Ionomer Composite Cement for lining under all types of material. Can also be used as an extended fissure sealant and for restoration of smaller lesions. Fluoride release and high biocompatibility properties to prevent secondary caries. High compressive strength (210MPa) and Radiopaque.

Ionoseal Glass Ionomer Composite Liner 2x4g tubes

VOC1126 

BASE LINER



**Amalgam Liner**  
VOCO

Is a ready to use varnish with silver to insulate dentine against chemical and bacterial influences. It is suitable for dentine sealing under amalgam, sealing of cements before applying amalgam and the sealing of amalgam fillings immediately after application.

Clinic Bottle 9g

VOC1026

CEMENTS



**Implant Cement**  
PREMIER

Premier Implant Cement offers the strength and durability of resin cements, ensuring long term performance for permanent implant retained restorations and long term provisionals.

Premier Implant Cement

Contains: 5ml Automix Dual Barrel Syringe, 10 Mixing Tips, 1 x Mixing Pad

PRE1450

Premier Implant Cement Refills

Contains: 25 Automix Dual Barrel Syringe Mixing Tips

PRE1451

FLOWABLE COMPOSITE



**Grandio® Flow**  
VOCO

A flowable universal light curing nano-hybrid restorative offering the same outstanding properties as Grandio.

Grandio Flow Refills – 2 x 2 g syringes

A1 VOC1863

A2 VOC1864

A3 VOC1865

A3.5 VOC1866

Also available in sets.



**Grandio Flow Caps**  
VOCO

Grandio Flow in caps for direct application.

Grandio Flow Caps Refills – 20 x 0.2s g caps

A1 VOC1880

A2 VOC1881

A3 VOC1882

A3.5 VOC1883

Also available in sets.



**Admira Flow**  
VOCO

A flowable light-curing Ormocer®-based restorative.

Admira Flow Refills – 2 x 1.8 g syringes

A1 VOC2486

A2 VOC2482

A3 VOC2483

A3.5 VOC2484



**Admira Flow Caps**  
VOCO

Admira Flow in Caps for direct application.

Admir Flow Caps – 25 x 0.25g, shades.

A1 VOC2491

A2 VOC2492

A3 VOC2493

A3.5 VOC2494

**FLOWABLE COMPOSITE**



**PermaFlo®**  
ULTRADENT



PermaFlo is a light cured, radiopaque, methacrylate-based flowable composite material. This flowable composite is highly thixotropic and exhibits very low film thickness compared with other flowable materials in the market. PermaFlo flows under pressure, but holds its shape and stays in place prior to light curing. Provides superb control of the material to prevent air entrapment while adapting the composite to cavity walls and undercuts. Ultradent Permaflo has superior high gloss polish and excellent properties to resist fracture and wear - ideal to use for both anterior and posterior restorations. The radiopaque property of PermaFlo will provide dental professionals trouble free future radiographic interpretations.

Universal Kit and Refill Kits contain:  
1 x 1.2ml of each Shade A1, A2, A3, A3.5, A4, B1, Dentin Opaquer, Translucent, 1 x 1.2 ml Syringe UltraEtch, 1 x 1.2ml Syringe PQ1 Bond, 1 x Half Size Syringe Organiser, 20 x Black Mini Tips, 20 x Micro 20gm Tips, 1 x Shade Guide.

ULT1273

Refills contain:  
2 x 1.2ml Syringe Packaging

Permaflo®		
SHADE CODE	BARCODE	
A1	ULT0947	
A2	ULT0948	
A3	ULT0949	
A3.5	ULT0952	
A4	ULT0954	
B1	ULT0956	
Translucent		
	ULT0612	
Dentin Opaquer		
	ULT1005	
Pink	ULT0963	
Purple	ULT0962	

**UNIVERSAL COMPOSITE**



**Grandio®**  
VOCO

A universal light cured nano-hybrid composite. The unique filler structure has a high filler degree (87%) that provides extreme hardness, very low shrinkage and excellent abrasion resistance. It can be use with all enamel and dentine bonding agents.

Grandio® – Syringe Refills		
A1	VOC1810	
A2	VOC1811	
A3	VOC1812	
A3.5	VOC1813	
B2	VOC1815	
C3	VOC1818	

Grandio® – Capsules		
A1	VOC1840	
A2	VOC1841	
A3	VOC1842	
A3.5	VOC1843	
C2	VOC1847	



UNIVERSAL COMPOSITE



Amelogen® Plus  
ULTRADENT

Amelogen Plus is a state-of-the-art Bis-GMA, radiopaque material for Class I, II, III, IV, V, VI restorations and direct veneers. It possesses wear resistance, strength and polishability making this simple composite system the perfect material for both posterior and anterior restorations. It also displays exceptional handling, optical and polishing characteristics. The unique syringe design of Amelogen Plus features the QuadraSpense and Kleen Sleeve. The QuadraSpense makes it possible to remove small amounts of material from the syringe. The Kleen Sleeve ensures that no black specks from the outer syringe become embedded in the composite.

Amelogen Plus Syringe Complete Kit:  
15 Shades- 1 x 2.5 grams syringe of A1,A2,A3, A3.5, A4, A5, B1, C2 Opaque White, Enamel White, Enamel Neutral, Enamel Gray, Translucent White, Translucent Gray, Translucent Orange 1- Shade guide, 1 full syringe organizer, 1 x 1.2 ml PQ1, 1 x 1.2 ml UltraEtch, 1 x 1.2ml PermaSeal

ULT3254

Amelogen Plus Singles Complete Kit:  
15 Shades- 1 x 10 pk of 0.3gm of each shade A1,A2,A3, A3.5, A4, A5, B1, C2 Opaque White, Enamel White, Enamel Neutral, Enamel Gray, Translucent White, Translucent Gray, Translucent Orange 1- Shade guide, 1 full syringe organizer, 1 x 1.2 ml PQ1, 1 x 1.2 ml UltraEtch, 1 x 1.2ml PermaSeal

ULT4525

Amelogen Plus Syringe Cosmetic Kit:  
7 Shades- 1 x 2.5g A1,A2,A3,Opaque White, Enamel White, Enamel Neutral, Translucent White, 1- Shade guide, 1 syringe organizer, 1.2 ml PQ1, UltraEtch, PermaSeal

ULT0315

Amelogen Plus Singles Cosmetic Kit:  
7 Shades- 1 x 10 pack singles A1,A2,A3,Opaque White, Enamel White, Enamel Neutral.

ULT3104

Amelogen Plus Syringe Basic Kit:  
7 Shades- 1 x 2.5g A1,A2,A3,A4,A5, B1, C2, 1- Shade guide, 1 syringe organizer, 1.2 ml PQ1, UltraEtch, PermaSeal.

ULT3098

Amelogen Plus Singles Basic Kit:  
7 Shades- 1 x 10 pack singles A1,A2,A3,A4,A5, B1, C2,Shade guide, 1 syringe organizer, 1.2 ml PQ1, UltraEtch, PermaSeal.

ULT3209

Amelogen® Plus Syringe Refills - Pack of 1 x 2.5gm syringe.

SHADE	CODE	BARCODE
A1	ULT9030	[Barcode]
A2	ULT9031	[Barcode]
A3	ULT9032	[Barcode]
A3.5	ULT9028	[Barcode]
A4	ULT9033	[Barcode]
A5	ULT9034	[Barcode]
B1	ULT9035	[Barcode]
C2	ULT9036	[Barcode]
Opaque White	ULT9037	[Barcode]
Enamel White	ULT9038	[Barcode]
Enamel Neutral	ULT9039	[Barcode]
Enamel Grey	ULT9040	[Barcode]
Translucent White	ULT9041	[Barcode]
Translucent Gray	ULT9042	[Barcode]
Translucent Orange	ULT9043	[Barcode]

Amelogen® Plus Singles Refills® - Pack of 10 x .3gm syringe.

SHADE	CODE	BARCODE
A1	ULT8010	[Barcode]
A2	ULT8011	[Barcode]
A3	ULT8013	[Barcode]
A3.5	ULT8024	[Barcode]
A4	ULT8014	[Barcode]
A5	ULT8015	[Barcode]
B1	ULT8017	[Barcode]
C2	ULT8020	[Barcode]
Opaque White	ULT8021	[Barcode]
Enamel White	ULT8022	[Barcode]
Enamel Neutral	ULT8023	[Barcode]
Enamel Gray	ULT8012	[Barcode]
Translucent White	ULT8016	[Barcode]
Translucent Gray	ULT8018	[Barcode]
Translucent Orange	ULT8019	[Barcode]

UNIVERSAL COMPOSITE



Vit~L~escence™  
ULTRADENT

The strength and wear resistance of a hybrid with the beauty and polishability of a microfill. Bis-GMA resin suitable for anterior and posterior restorations, Vit~L~escence offers superior fluorescence and opalescence allowing restorations to mimic natural enamel. With low polymerisation shrinkage and high abrasion resistance and colour stability, Vit~L~escence offers a new and unique concept for exquisite restorations. 2.5gm syringes.

Vit~L~escence Masters Kit:

Contains: 1 x syringe of all Vit~L~escence shades, complete shade guide, full size syringe organiser














ULT0500 

Vit~L~escence Advance Kit:

Contains: 1 x syringe of A1, A2, A3, A4, A5, B1, OS, PF, PN, PA, PS, TM, IB, Advance shade guide, full size syringe organiser

ULT0515 

Vit~L~escence™ - Dentin Refill

SHADE	CODE	BARCODE
A1	ULT0358	
A2	ULT0343	
A3	ULT0344	
A3.5	ULT0356	
A4	ULT0360	
A5	ULT0362	
A6	ULT0408	
B1	ULT0409	
B2	ULT0418	
B3	ULT0421	
B4	ULT0422	
B5	ULT0423	
C1	ULT0426	
C2	ULT0435	
C3	ULT0439	
C4	ULT0440	
C5	ULT0441	
D3	ULT0442	

Vit~L~escence™ Enamel Refill

SHADE	CODE	BARCODE
Translucent Grey	ULT0478	
Translucent Ice	ULT0479	
Translucent Mist	ULT0482	
Translucent Smoke	ULT0485	
Translucent Yellow	ULT0486	
Translucent Amber	ULT0499	
Translucent Frost	ULT1187	
Translucent Blue	ULT1188	
Translucent Orange	ULT1189	
Opaque White	ULT1182	
Opaque Snow	ULT1183	
Pearl Neutral	ULT1184	
Pearl Amber	ULT1185	
Pearl Smoke	ULT1186	
Iridescent Blue	ULT1317	
Pearl Frost	ULT0443	
Vitalescence Masters Shade Guide	ULT0822	
Vitalescence Advanced Shade Guide	ULT1171	

UNIVERSAL COMPOSITE



**THE DENTAL ADVISOR**  
++++½

**Amaris®**  
VOCO

Highly aesthetic light-curing composite.

Amaris Syringe Set: (8x4g syringes 01, 02, 03, 04, 05, TL, TN, TD/ 2x 1.8g syringes Amaris Flow-HT, HO) & Accessories  
VOC1930

Amaris Caps Set: (16x0.25 g caps-01, 02, 03, 04, 05, TL, TN, TD/ 2x1.8g syringes Amaris Flow(HT,HO) & Accessories  
VOC1950

Amaris Opaque Syringe Refills  
1 x 4g syringe in various shades

01	VOC1933	
02	VOC1934	
03	VOC1935	
04	VOC1936	
05	VOC1937	
0	VOC1938	
Translucent		
TL	VOC1940	
TN	VOC1941	
TD	VOC1942	

Amaris Opaque Caps Refills  
16x 0.25 g Caps in various shades

01	VOC1953	
02	VOC1954	
03	VOC1955	
04	VOC1956	
05	VOC1957	
0	VOC1958	
Translucent		
TL	VOC1960	
TN	VOC1961	
TD	VOC1962	



**THE DENTAL ADVISOR**  
++++  
**Arabesk Top / Caps**  
VOCO

Light-curing microhybrid composite with BCS (Bi-Ceramic System) filler system.

Arabesk Set: (7x 4g syringes A2, A3, A3.5, B2, B3, C2, OA2), 5ml Vococid gel, 4ml solobond M, accessories.

VOC1530

Arabesk Syringe Refill - 4g

A1	VOC1532	
A2	VOC1533	
A3	VOC1534	
A3.5	VOC1535	

POSTERIOR COMPOSITE



**Clinical Case**



Source:  
Dr. Ulrich Sachse, Kiel

**X-tra fil**  
VOCO

Light-cured posterior filling material for restoration of class I and II in bulk filling technique with a curing depth of 4mm and a curing time of only 10s. Can also be used as a core build-up. Available in caps for direct application.

X-tra fil Composite Universal Capsules 20x 0.25gm  
VOC1741

X-tra fil Composite Universal Syringe 10 x 5gms  
VOC1739

X-tra fil Composite Universal Syringe x 5gms  
VOC1740

Easygloss  
VOC2050

TEMPORARY FILLING



Clip  
VOCO

Light-curing temporary filling material. Suitable for temporary fillings, temporary sealing of cavities, especially in inlay/onlay technique.

- Easy to place - easy to remove in one piece
- No damage to preparation boundaries
- Time-saving and economical
- Tight margins and durable

Syringes 2x 4 g

VOC1284 

Tripack Syringes 3x 4 g

VOC1285 



Clip F  
VOCO

A fluoride releasing Clip - fluorides are effective against secondary caries.

Syringes 2x 4 g

VOC1289 

Tripack Syringes 3x 4 g

VOC1283 



Clip Caps  
VOCO

Clip in practical Caps for direct time-saving application.

Caps 25x 0.25 g

VOC1287 

Caps Dispenser, 1 piece

VOC9301 

GLASS IONOMER



Twinky Star  
VOCO

A coloured light-cured compomer with glitter effect for the restoration of deciduous teeth. Easily and quickly applied from capsules, it has outstanding polishability and fluoride release.

Kit contains:

- 40 x 0.25g capsules in different colours, & shade guide

VOC1680 

Refills 25 x 0.25g, gold

VOC1681 

Refills 25 x 0.25g, blue

VOC1683 

Refills 25 x 0.25g, green

VOC1685 

Refills 25 x 0.25g, lemon

VOC1687 



Meron Plus AC  
VOCO

Universal resin-reinforced glass ionomer luting cement in application capsules. Self-adhesive: saves time - no separate primer or conditioner needed. Thixotropic: no unwanted run-off. Specifically suited for zirconia-based pieces.

Pack of 48 capsules

VOC1733 



Ionofil Molar AC  
VOCO

A packable glass ionomer cement in application capsules for direct application into molar fillings and into core build ups. It has 40% more available cement in each capsule compared with similar products and provides safe adhesion to dentine and enamel.

Pack 48 capsules

A1 VOC1464 

A3 VOC1465 

B3 VOC1466 

VOC1466 



Ionofil Molar AC Quick  
VOCO

Is an immediately packable glass ionomer cement in application capsules. Suitable for long term temporary class I and II restorations. It has a setting time of only 2.5 minutes making it the perfect material for the treatment of children.

Starter Kit contains:

- 48 capsules (8 x A1, 32 x A3, 8 x B3) & 3 ml Final Varnish LC

VOC1630 

Refills Packs 48

A1 VOC1631 

A2 VOC1634 

A3 VOC1632 

VOC1632 

CORE BUILD UP MATERIAL



**CompCore AF**  
PREMIER

CompCore AF is a composite specifically designed for Core build-ups. Its patented hyperbranched molecule promotes a dense cross-linking of composite monomers for more complete curing, added strength and reduced shrinkage.



**Refill Kits**  
PREMIER

Refill Kits Contain:  
50 gm cartridge  
25 x syringes  
25 x mixing tips  
25 x syringes  
25 x mixing tips & 10 x Accor matrices

White  
PRE1432

A3  
PRE1433

Applicator Gun Type-25  
PRE1416

Mixing Tips and Intraoral Tips (35 each)  
PRE1417



**CompCore AF Stack**  
PREMIER

AF Stack SyringeMix White  
PRE1512

AF Stack SyringeMix A3  
PRE1513



**Rebilda DC**  
VOCO



A dual cure flowable core build up system for use on both vital and non-vital teeth. Low heat generation offers safety during bulk applications while the optimised small mixing tips saves on material.

Dentine Kit Contains: 50g cartridge dentine, 4ml each Futurabond DC Liquid 1 and 2, mixing tips type 12, dispenser type 3, intraoral tips type 3, accessories

VOC1395

Refill Kits Contain: refill cartridge 50g, Rebilda DC dentine, mixing tips type 12, intraoral tips type 3

VOC1396

Refill Kits with Rebilda DC blue  
VOC1397

Refill Kits with Rebilda DC white  
VOC1398



**Rebilda DC - QuickMix Set**

QuickMix Set: 10g syringe dentin, 4ml each Futurebond DC liquid 1 & 2, mixing tips type 11, intraoral tips type 4, accessories.

Quickmix Syringe Set  
VOC1402

QuickMix Refill Dentine  
VOC1403

Quickmix Refill Blue  
VOC1404

Quickmix Refill White  
VOC1405

MATRICES & STRIPS



**Omni-Matrix™**  
ULTRADENT

Omni-Matrix with wingless design provides placement of wedges with ease. Bands are colour coded for easy and fast selection. Omni-Matrix band is thin, adaptive and burnishable. Adjust the band's circumference by twisting the conical handle. Its patented, articulated head swivels, fitting comfortably into any quadrant of the mouth. Once restoration is completed, the Omni-Matrix is easily released, leaving restorative material intact. Equipped with an internal screw mechanism attached to a conical handle that adjusts the band's circumference to adapt to the tooth. The patented rotating head means your bonding is always in correct orientation.

Intro Kit Contains: 16 x 0.001 and mylar bands  
8 x 0.0015 and pedo bands

ULT1105

**Refills**

0.001 (Orange) – Pack 48

ULT1101

0.0015 (Green) – Pack 48

ULT1102

Mylar (Red) – Pack 48

ULT1103

Pedo (Purple) – Pack 48

ULT1104



**Waterpik Tofflemire Bands**  
PREVIDENT

No.1 Adult Universal (0.002) – Pack 36

PVDWP345770

No.1 Adult Universal (0.0015) – Pack 36

PVDWP345789

No.2 Adult Universal (0.001) – Pack 36

PVDWP346044

No.2 Adult Mod Wide(0.002) – Pack 36

PVDWP345797

No.2 Adult Mod Wide (0.0015) – Pack 36

PVDWP345800

No.1 Adult Universal Dead Soft (0.001) – Pack 36

PVDWP346036

No.13 Pedo Universal (.001 – Pack 36)

PVDWP346052

No.13 Pedo Universal (.002)

PVDWP345842

No.13 Pedo Universal (.0015)

PVDWP345850



**Matrix Bands & Dispensers**

FRASACO

Plastic universal strip for use with self cure and light cure materials. They are available as straight strips, or rolls which fit into a patented matrix band dispenser. Available in transparent, or yellowish transparent, although not recommended for light curing, it makes the band easier to see on the tooth or the cabinet top.

Universal Strips: Straight, L: 120mm, W: 8.5mm, 50 per pack. Transparent, extra-thin (0.05mm)  
FRAUS-120KE

Matrix Rolls: Matrix band dispenser with one roll L: 15m, W: 8.5mm. Transparent, extra thin: 0.05mm. Pkt 3 Rolls.  
FRAUS-15PKE

Matrix Band Refill Roll: L: 15mm, W: 8.5mm, 3 rolls per pack. Transparent extra-thin, 0.05mm  
FRAUS-15PKEN

Yellowish-Transparent extra-thin, 0.05mm  
FRAUS-15PEEN



**InterGuard®**  
ULTRADENT

InterGuard ensures a faster, safer preparation by protecting the adjacent tooth from iatrogenic damage. Durable stainless steel InterGuard is autoclavable and can be used multiple times. Stable curls at each end leave transition angles clear for full access. Great for tunnel preparations and protecting the adjacent tooth during air abrasion.

Refills: InterGuard 4.0mm. 10pk.

ULT4011

Refills: InterGuard 5.5mm. 10pk.

ULT4012



**Quantum Separating Strips**

ADM

Quantum mylar separating strips in roll form provide the balance of quality and economy. Available in ultrathin mylar (0.002' thick) in 8mm and 10mm widths, 15 metres in length along with a stainless steel cutter integrated within the dispensing box. They suit all requirements and ensure tight contacts on composites.

6mm wide, Clear

QTM06

8mm wide, Clear

QTM08

10mm wide, Clear

QTM10

**MATRICES & STRIPS**



**Tofflemire Bands**

PD

Stainless steel, thickness: thin - 0.05mm, extra thin: 0.035mm. Packs of 12.

- Thin 01
- PDTBT01
- Thin 02
- PDTBT02
- Thin 03
- PDTBT03
- Thin 11
- PDTBT11
- Thin 13
- PDTBT13
- Thin Assorted
- PDTBTASS
- X-Thin 01
- PDTBXT01
- X-Thin 02
- PDTBXT02
- X-Thin 03
- PDTBXT03
- X-Thin 11
- PDTBXT11
- X-Thin 13
- PDTBXT13
- X-Thin 14
- PDTBXT14
- X-Thin 15
- PDTBXT15



- Junior Retainer – With 12 Bands
- PDTRJ010
- Adult Retainer – With 12 Bands
- PDTRS011
- Senior Retainer – With 12 Bands
- PDTBXT13
- Universal Retainer – Pack 12
- PDTRU012

**WEDGES**



**Interdental**

PROGRESS

- Extra Thin Neutral – Pack 250
- PRO920
- Thin Neutral – Pack 250
- PRO921



**Wizard Wedges**

PREIDENT

Finest quality white birchwood; will not slip or splinter. Packs of 500.

- Small
- PVD61211
- Large
- PVD61212
- Slim Jim
- PVD61213
- Wide Base
- PVD61214
- Assorted
- PVD61215

**ARTICULATING PRODUCTS**



**Paper**

PD

- Thin Blue - Pack 12 (50 Microns)
- PDAPB001
- Thick Blue - Pack 12 (80 microns)
- PDAPB002



**Booklet**

BAUSCH

At 40 microns these Articulating papers are extremely thin and tear resistant and are coated with liquid colours on both sides. They are especially suitable for the two tone representation of static and dynamic occlusion.

- Blue
- DIMBK09
- Red
- DIMBK10



**Dispenser**

BAUSCH

- Plastic Dispenser Blue
- DIMBK01
- Plastic Dispenser Horseshoe Blue
- DIMBK03

ARTICULATING PRODUCTS



**Progress 100**  
BAUSCH

This new smooth fiber reinforced paper adapts perfectly to the occlusal surfaces. The contact points are marked very precisely because of its progressive colour transfer. 100 micron thin paper which is impregnated with hydrophile waxes and pharmaceutical oils.

Plastic Dispenser 300 Strips Blue  
DIMBK51

Plastic Dispenser 300 Strips Red  
DIMBK52



**Foil**  
BAUSCH

Ultra Thin 8 microns. Especially suited for representing static and dynamic occlusion.

Arti-Fol - 22mm 20m Roll  
DIMBK27



**Forcep**  
BAUSCH


Miller Style  
DIMBK142



**Articulating Forceps**  
PREMIER

PRE65130

FINISHING & POLISHING



**Jiffy Polishes®**  
ULTRADENT

Abrasive, latex-free material impregnated silicon rubber cups, disks and points of appropriate grits and shapes excellent for contouring, finishing, and polishing composites such as Amelogen Plus or Vit-I-essence. Packs of 20.

Jiffy Polisher Cups Coarse (green)  
ULT0890

Jiffy Polisher Disks Coarse (green)  
ULT0891

Jiffy Polisher Points Coarse (green)  
ULT0892

Jiffy Polisher Cups Medium (yellow)  
ULT0838

Jiffy Polisher Points Medium (yellow)  
ULT0839

Jiffy Bur Polisher Disks Medium (yellow)  
ULT0840

Jiffy Polisher Cups Fine (White)  
ULT0841

Jiffy Polisher Points Fine (White)  
ULT0842

Jiffy Polisher Disks Fine (White)  
ULT0843



**Unique Polishes**  
NTI

Slow speed Single Step composite polishing system for high lustre polishing of all types of composite resin. The new silicone / diamond impregnation technology of Unique Polishes means that the perfect grit is always selected for the composite being polished. The stainless steel gold coloured shank provides true concentricity and strength as compared with other plastic type shanks and chairside polishing time is greatly reduced.

Point. Available in Packs of 6.

NTI-P20032

Cups in Packs of 6.

NTI-P20035

Disc in Packs of 6.

NTI-P20038



**Diamond Twist SCL & SCO**  
PREMIER

Create a glaze like finish on castable ceramics and CAD/CAM type porcelain or composite crowns. Diamond Twist's high viscosity, densely concentrated diamond polish paste and has 130% more diamonds with optimum sized crystals for dramatically improved ceramic polishing results.

Diamond Polish SCL – Extra-Oral  
PRE19000

Diamond Polish SCO – Intra-Oral  
PRE19002

FINISHING & POLISHING



**Compo-Disc**  
PREMIER

Hand held diamond Finishing Disc. Use Compo-Disc to open proximal contacts and for contouring. Thin, Single-sided. Edge and face cutting.

Compo-Disc (pkg of 2 discs)  
PRE13106



100T  
PRE13107  
100F  
PRE13111  
100UF  
PRE13115  
150T  
PRE13109  
150F  
PRE13113  
150UF  
PRE13117

Grit	Thickness
T = 60 $\mu$	(0.15mm)
F = 45 $\mu$	(0.127mm)
UF = 20 $\mu$	(0.1mm)

**Compo-Strips®**  
PREMIER

These thin flexible one-sided diamond abrasive strips have an uncoated safe centre for easy access through the contact point and are ideal for finishing class III, IV and veneer restorations.

- 100T – Blue. Pack 6  
PRE13107
- 150T – Blue. Pack 6  
PRE13109
- 100F – Red. Pack 6  
PRE13111
- 150F – Red. Pack 6  
PRE13113
- 100UF – Yellow. Pack 6  
PRE13115
- 150UF – Yellow. Pack 6  
PRE13117
- Assorted Strips – Pack 6  
PRE13127



**Serrated Diamond Finishing Strips**  
NTI

Single sided – safe for approximal teeth. Quality raw diamonds provide smooth and long-lasting use. Serrated mid-section to help tackle tight contact points. Autoclavable, available in 3 grits – Medium, Fine and Extra Fine. Metal thickness ranges from 0.13 – 0.08mm.

- Medium – 37-44 microns. Pack of 10  
NTI-FS3MX
- Fine – 30-40 microns. Pack of 10  
NTI-FS3FX
- Superfine – 0-20 microns. Pack of 10  
NTI-FS3SFX



- Wide Coated on one side. Pack of 10  
NTI-FS3X
- Wide Coated on one side. Pack of 10  
NTI-FS3MX



**Luminescence™ Plus**  
PREMIER

Direct and indirect diamond polishing paste that contains Potassium Nitrate for patient comfort.

Luminescence Plus Gel (3gm)  
PRE14200



**Diamond Polish Mint**  
ULTRADENT

Fine Microcrystalline diamond particles for unsurpassed polishing of all porcelain and composite restorations. Mint flavour for patient satisfaction.

Diamond Polish Refill - 2 x 1.2ml syringes  
ULT5541

COMPOSITE WETTING AGENTS



**Composite Wetting Resin**  
ULTRADENT

Composite Wetting Resin is a 45% filled, light cured, liquid resin. Significantly superior to the use of single-component adhesives, which contain solvents and inhibit composite polymerisation. Moistens dry composite during contouring.

Composite Wetting Resin - 2 x 1.2ml syringes  
ULT3059

COMPOSITE SEALER & VARNISH

DESENSITISER



**PermaSeal®**  
ULTRADENT

A composite sealer and bonding agent that is a light cured, methacrylate based, unfilled resin. It cures with a very low film thickness and reduces microleakage when placed on class V composite margins. It should be applied immediately upon completion of a restoration so it can penetrate and seal, thus reducing wear and marginal breakdown.

Kit - 4 x 1.2ml syringes & 20 Tips

ULT0630

Refill - 4 x 1.2ml syringes

ULT0631



**Enamel Pro® Varnish Clear**  
PREMIER

The ONLY fluoride varnish formulated to deliver ACP (Amorphous Calcium Phosphate). ACP works together with fluoride to desensitise the dentine by providing a semi-permanent occlusion with fluorapatite. Enamel Pro Varnish CLEAR with ACP technology has been proven to deliver 4X more fluoride uptake to enamel than any other transparent varnish on the market.

Each box contains single-unit doses and an applicator brush, one Comfy-grip™ cup holder, Post-application patient care instructions and doctor's directions for use.

Clear Strawberry'n'Cream 0.40 ml Pkt 35

PRE7540

Clear Strawberry'n'Cream 0.40 ml Pkt 200

PRE7543

Clear Bubblegum 0.40ml Pkt 35

PRE7541

Clear BubbleGum 0.40ml Pkt 200

PRE7544

Clear BubbleGum 0.25ml Pkt 35

PRE7542



**Bifluorid 10**  
VOCO

Fluoride varnish for the treatment of hypersensitivity.

Bifluorid 10 1x 10gram Bottle.

VOC1616



**Admira® Protect**  
VOCO

A light cured Ormocer based desensitiser with special filler technology for high abrasion proofness. Suitable for the treatment of hypersensitive dentine in cervical areas and crown margins.

Pack 4 x 5ml bottles & accessories

VOC1650

SingleDose pack of 50

VOC1651

SingleDose pack of 200

VOC1652

Brands you can trust

- ★ BAUSCH
- ★ NTI
- ★ PD
- ★ PIERRE ROLLAND
- ★ PREMIER
- ★ PREVIDENT
- ★ RIBBOND
- ★ SATELEC
- ★ ULTRADENT
- ★ VOCO

For all your restorative needs

TALK TO US

Opalescence<sup>®</sup>  
tooth whitening systems



Be Dazzling



Opalescence offers whitening options for every lifestyle and budget. With a PF formula that helps to improve overall enamel health, and a unique, viscous gel that safely and effectively whitens teeth, Opalescence creates dazzling smiles.

Be Dazzling with Opalescence.

 **ULTRADENT**  
PRODUCTS, INC.

phone: 0800 30 10 10

fax: 0800 486 9336

[www.gunz.co.nz](http://www.gunz.co.nz)