

Diamantschleifinstrumente • Diamond Grinding Instruments • Fraises diamantées de laboratoire

Labor • Laboratory • Laboratoire

379
(277)

Eiform
Egg
Oeuf



Größe/Size/Taille $\varnothing$ 1/10 mm	014
Länge-Length in mm - Longueur	3,0
Shank   	
HP  806 104 277 534 ... 379 - ... C - HP	014
HP  806 104 277 514 ... 379 - ... F - HP	014
 opt. 30.000 rpm	

379
(277)

Eiform
Egg
Oeuf



Größe/Size/Taille $\varnothing$ 1/10 mm	023
Länge-Length in mm - Longueur	4,2
Shank   	
HP 806 104 277 524 ... 379 - ... M - HPA	023
 opt. 20.000 rpm	

801
(001)





Rund
Round
Rond



Größe/Size/Taille $\varnothing$ 1/10 mm	009	010	012	014	016	018	021	023	029	035
Länge-Length in mm - Longueur										
Shank   										
HP 806 104 001 524 ... 801 - ... M - HP	009	010	012	014	016	018	021	023	029	035
 opt. 30.000 rpm, > 023-035 20.000 rpm										

805
(010)

Umgekehrter Kegel
Inverted Cone
Cône renversé



Größe/Size/Taille $\varnothing$ 1/10 mm	012	014	016	018	021	023	040	050
Länge-Length in mm - Longueur	1,5	1,5	1,5	1,5	2,0	2,0	1,5	2,0
Shank   								
HP 806 104 010 524 ... 805 - ... M - HP	012	014	016	018	021	023	040	050
 opt. 40.000 rpm, > 023-035 30.000 rpm, > 035 20.000 rpm								

807
(225)

Umgekehrter Kegel
Inverted Cone
Cône renversé



Größe/Size/Taille $\varnothing$ 1/10 mm	016	018	023
Länge-Length in mm - Longueur	4,0	5,0	6,0
Shank   			
HP 806 104 225 524 ... 807 - ... M - HP	016	018	023
 opt. 30.000 rpm			



805-014M-HP

808
(014)

Umgekehrter Kegel
Inverted Cone
Cône renversé



Größe/Size/Taille $\varnothing$ 1/10 mm	018	023
Länge-Length in mm - Longueur	0,7	0,8
Shank   		
HP 806 104 014 524 ... 808 - ... M - HP	018	023
 opt. 30.000 rpm		

809
(019)

Umgekehrter Kegel mit Ansatz
Inverted Cone with Collar
Cône renversé avec extension



Größe/Size/Taille $\varnothing$ 1/10 mm	023
Länge-Length in mm - Longueur	3,5
Shank   	
HP 806 104 019 524 ... 809 - ... M - HP	023
 opt. 30.000 rpm	

Bitte Drehzahlbereiche beachten (Drehzahltable - Seite 11.14) / Please observe Recommended speeds (speed chart - page 11.14) / Respecter les vitesses de rotation recommandées (voir à la page 11.14)

 F	Roter Ring/Red ring/Bague rouge	Fein/Fine/Fin	50 μ m	ISO No. 514
 M	Ohne Ring/Without ring/Sans bague	Standard/Standard/Standard	110 -120 μ m	ISO No. 524
 C	Grüner Ring/Green ring/Bague verte	Grob/Coarse/Gros	135 μ m	ISO No. 534
 SC	Schwarzer Ring/Black ring/Bague noire	Supergrob/Supercourse/Supergros	180 μ m	ISO No. 544